# AKUS®

## SET ANESTESIA EPIDURAL UN SOLO USO

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ES

## EPIDURAL ANAESTHESIA SET FOR SINGLE USE

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#### EPIDURAL SET - INSTRUCTIONS FOR USE

##### Indications:

Localisation in the epidural space for a single puncture. Localisation and puncture of the epidural space, interpleural space, and interfascial spaces (for example, TAP block). In addition to transverse abdominal plane.  
The use of these types of needles is restricted to veterinary personnel with sufficient training and experience required to perform anaesthesia and/or epidural analgesia techniques as well as all other types of space or blocks.

##### Instructions for use:

1. Inform the patient owner about the technique and obtain written consent before starting the procedure.
2. **Closely** inspect the proposed puncture site for infection and lesions before performing the puncture.
3. Position the animal properly, according to the technique to be carried out. Correctly disinfect the area and prepare a sterile drape before starting the procedure. To ensure the accuracy of the block and animal welfare, it is recommended to perform the procedure after pre-medication or prior sedation.
4. Inspect the puncture material **closely**, ensuring that it has been packaged correctly and the needle is sterile. If not in perfect condition or sterility has been compromised, discard immediately.
5. The epidural needle must be inserted with the catheter placed inside. It may be used with or without wings. The large, ergonomic AKUS™ wings allow for greater accuracy in the insertion of the metal cannula. Once it is 2-3 cm inside, guided by the centimetre markings on the cannula, remove the stylet and insert the low-resistance syringe, which can be used with air or 0.9% saline solution for the loss-of-resistance technique.
6. Insert the needle **slowly** until the epidural space is located by the loss of resistance with air or serum. Check the correct location by means of the loss of resistance technique, hanging drop or Baraka test.
7. Once the needle is in the correct position, insert the medication.

##### The set includes:

- Tuohy epidural needle tip.
- Low-resistance syringe to locate the epidural space using the LOR technique.
- AKUS SAFETY device.

##### Contraindications:

1. **Absolute:**
  - Serious coagulation disorders, anti-coagulant treatment.
  - Skin infection in the place of puncture.
  - Septicaemia.
  - Allergy to local anaesthetics.
  - Severe intracranial hypertension.
2. **Relative:**
  - Anatomical malformation.
  - Skin infections distant from the puncture site.
  - Mild coagulations disorders.
  - Uncorrected blood volume deficit.
  - Dehydration.

##### Procedure complications:

- Inadvertent puncture of the intrathecal space.
- PDPH (Postdural puncture headache) – post lumbar puncture headache.
- Horner's Syndrome.
- Hypotension.
- Urinary retention.
- Lower back pain.
- Paresthesia.
- Partial paralysis.
- Temporal central block.
- Subdural or epidural haematoma.
- Meningitis.
- Systemic toxic reaction.
- Central block.
- Infection (epidural abscess).
- Respiratory depression.
- Cardiac arrest.
- Death.

##### Warning:

- **All** AKUS™ epidural kits have been sterilised with ethylene oxide.
- Single-use needle and material. Do not reuse.
- The AKUS™ needle is phthalate-free.
- Do not use if the packaging has been damaged.
- Ensure that the necessary material for general anaesthesia and cardiopulmonary resuscitation is readily available.